

SIEMENS DIMENSION EXL, RXL & Xpand SPECIFIC PACKAGING

**UNIQUE ASSAY AND PACKAGING SOLUTIONS
TO FIT YOUR LABORATORY NEEDS!**



All Products 510(k) Cleared

THREE VIAL ENZYMATIC HOMOCYSTEINE

- No significant interference from cystathionine which is present in renal disease
- Excellent correlation to HPLC and immunochemical methods
- Validated application parameters available for the Xpand

DIRECT ENZYMATIC HbA1c

- Diazyme's Enzymatic HbA1c virtually eliminates interference from hemoglobin variants
- The assay's unique enzymatic methodology directly measures glycosylated hemoglobin
- The assay correlates with Tosoh HPLC and Roche Tina-Quant methods
- Eliminates latex particle build up, thereby reducing instrument maintenance requirements
- Validated application parameters available for the Xpand
- Suggested application parameters available for EXL

CYSTATIN C

- Calibrator Traceability (ERM-DA471/IFCC)
- Improved on-board and calibration stability
- Validated application parameters available for the Xpand
- Suggested application parameters available for EXL

