

DZ-LITE™ c270 INSTRUMENT

Assay Menu Includes

IVD


- Apolipoprotein A-1
- Apolipoprotein B
- Beta-Hydroxybutyrate◆
- Carbon Dioxide Enzymatic
- Creatinine Liquid Reagent
- Cystatin C
- D-Dimer
- Direct Enzymatic HbA1c
- Direct HbA1c
(Enzymatic, On-Board Lysis)
- Enzymatic Lithium
- Enzymatic Potassium
- Enzymatic Sodium

IVD (continued)

- EZ Vitamin D
- Glycated Serum Protein
- HDL-Cholesterol
- Homocysteine
- hsCRP
- Human Kappa FLC
- Human Lambda FLC
- LDL-Cholesterol
- Lipase◆
- Lipoprotein (a)
- PLAC® Test
for Lp-PLA2 Activity
- Procalcitonin
- Retinol Binding Protein◆
- Total Bile Acids
- 1,5-Anhydroglucitol (1,5-AG)
- 5'-Nucleotidase◆

USA: For Research Use Only*

- Adenosine Deaminase
- Adiponectin
- Cardiac Troponin I
- HDL3
- Heart-Type Fatty Acid-Binding Protein (H-FABP)
- Myeloperoxidase (MPO)
- N-Acetyl-β-Glucosaminidase (NAG)
- sdLDL



◆ Not CLIA Categorized

* USA: For Research Use Only. Not For Use In Diagnostic Procedures



PRODUCT OFFERINGS



Diazyme Europe GmbH

Zum Windkanal 21
01109 Dresden,
Germany
Tel: +49-351-886-3300
Fax: +49-351-886-3366
www.diazyme.de
sales@diazyme.de

Diazyme Laboratories, Inc.

12889 Gregg Court
Poway, CA 92064 USA
Tel: +1-858-455-4768
+1-888-DIAZYME
Fax: +1-858-455-3701
www.diazyme.com
sales@diazyme.com

Diazyme Shanghai Co., Ltd.

Room 201, 1011 Halei Road
Zhangjiang Hi-tech Park
Shanghai, 201203
People's Republic of China
Tel: +86-21-51320668
Fax: +86-21-51320663
www.lanyuanbio.com
service@lanyuanbio.com



Innovations In Clinical Diagnostics

CAT. #	CANCER MARKERS	METHOD	REGULATORY
DZ169A	Human Kappa Free Light Chain Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared Health Canada Registered
DZ170A	Human Lambda Free Light Chain Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared Health Canada Registered
DZ863A	Carcinoembryonic Antigen Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ USA: For Research Use Only*
DZ869A	Alpha-Fetoprotein Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ USA: For Research Use Only*
CAT. #	CARDIOVASCULAR MARKERS	METHOD	REGULATORY
DZ128A	LDL-Cholesterol Assay	(Dual Vial Liquid Stable) (PVS\PEGME Coupled)	☑️ 510(k) Cleared
DZ129A	HDL-Cholesterol Assay	(Dual Vial Liquid Stable) (PVS\PEGME Coupled)	☑️ 510(k) Cleared
DZ131B	Lipoprotein (a) Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared
DZ141A	Apolipoprotein A-1 Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared
DZ142A	Apolipoprotein B Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared
DZ145A	Cardiac Troponin I Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ USA: For Research Use Only*
DZ147A	Myoglobin Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared
DZ560A	Small Dense Low-Density Lipoprotein-Cholesterol (sdLDL) Assay	(Dual Vial Liquid Stable) (Enzymatic)	USA: For Research Use Only*
DZ561A	High-Density Lipoprotein III (HDL3) Assay	(Dual Vial Liquid Stable) (Enzymatic)	USA: For Research Use Only*
DZ565A	DZ-Lite OxPL-apoB Assay	(Immunochemiluminescent)	USA: For Research Use Only*
DZ568A	Homocysteine Enzymatic Assay	(Three Vial Liquid Stable) (Enzyme Cycling)	☑️ 510(k) Cleared Health Canada Registered
DZ568B	Homocysteine 2 Reagent Enzymatic Assay	(Dual Vial Liquid Stable) (Enzyme Cycling)	☑️ 510(k) Cleared Health Canada Registered
DZ785A	Heart-Type Fatty Acid-Binding Protein (H-FABP) Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ USA: For Research Use Only*
10-0148	PLAC® Test for Lp-PLA ₂ Activity	(Dual Vial Liquid Stable) (Enzymatic)	☑️ 510(k) Cleared
130206009M	hs-Cardiac Troponin I (hs-cTnI) Assay	(Immunochemiluminescent)	USA: For Research Use Only*
CAT. #	COAGULATION MARKERS	METHOD	REGULATORY
DZ179A	D-Dimer Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared
DZ768A	Fibrinogen Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	510(k) Cleared
CAT. #	DIABETIC MARKERS	METHOD	REGULATORY
DZ112B	Glycated Serum Protein Assay	(Dual Vial Liquid Stable) (Enzymatic)	☑️ 510(k) Cleared Health Canada Registered
DZ152A	1,5-Anhydroglucitol (1,5-AG) Assay	(Dual Vial Liquid Stable) (Enzymatic)	510(k) Cleared
DZ168A	Direct Enzymatic HbA1c Assay	(Liquid Stable) (Single Channel)	☑️ 510(k) Cleared Health Canada Registered
DZ168C	Direct HbA1c Assay (Enzymatic, On-Board Lysis)	(Dual Vial Liquid Stable) (Enzymatic)	☑️ 510(k) Cleared Health Canada Registered
DZ188A	Adiponectin Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	USA: For Research Use Only*
CAT. #	OTHER	METHOD	REGULATORY
DZ153A	Beta-Hydroxybutyrate Assay	(Dual Vial Liquid Stable) (Enzymatic Colorimetric)	510(k) Exempt
DZ187A	Retinol Binding Protein Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Exempt

CAT. #	ELECTROLYTE MARKERS	METHOD	REGULATORY
DZ113C	Enzymatic Potassium Assay	(Dual Vial Liquid Stable) (Enzymatic)	☑️ 510(k) Cleared
DZ114B	Enzymatic Sodium Assay	(Dual Vial Liquid Stable) (Enzymatic)	☑️ 510(k) Cleared
DZ116B	Enzymatic Lithium Assay	(Dual Vial Liquid Stable) (Enzymatic)	☑️ 510(k) Cleared Health Canada Registered
DZ122A	Carbon Dioxide Enzymatic Assay	(Single Vial Liquid Stable) (Enzyme Cycling 405nm)	☑️ 510(k) Cleared
CAT. #	INFECTIOUS DISEASE MARKERS	METHOD	REGULATORY
DZ901A	SARS-CoV-2 Neutralizing Antibody CLIA Kit	(Immunochemiluminescent)	☑️ USA: For Research Use Only*
130219015M	DZ-Lite SARS-CoV-2 IgG CLIA Kit	(Immunochemiluminescent)	FDA EUA
130219016M	DZ-Lite SARS-CoV-2 IgM CLIA Kit	(Immunochemiluminescent)	FDA EUA
CAT. #	INFLAMMATORY MARKERS	METHOD	REGULATORY
DZ135A	hsCRP Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared
DZ178C	Myeloperoxidase (MPO) Immunoassay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ USA: For Research Use Only*
DZ202A	Calprotectin Assay ♦	(Immunoturbidimetric)	USA: For Research Use Only*
CAT. #	LIVER MARKERS	METHOD	REGULATORY
DZ042A	Total Bile Acids Assay	(Dual Vial Liquid Stable) (Enzyme Cycling)	☑️ 510(k) Cleared Health Canada Registered
DZ082B	α-L-Fucosidase (AFU) Assay*	(Single Vial Liquid Stable) (Enzymatic)	☑️ USA: For Research Use Only*
DZ117A	Adenosine Deaminase Assay	(Dual Vial Liquid Stable) (Enzymatic)	☑️ USA: For Research Use Only*
DZ123A	5'-Nucleotidase Assay	(Dual Vial Liquid Stable) (Enzymatic)	☑️ 510(k) Exempt Health Canada Registered
CAT. #	LUNG DISEASE MARKERS	METHOD	REGULATORY
DZ902A	DZ-Lite MMP-1 CLIA Kit ♦	(Immunochemiluminescent)	USA: For Research Use Only*
DZ903A	DZ-Lite MMP-7 CLIA Kit ♦	(Immunochemiluminescent)	USA: For Research Use Only*
CAT. #	RENAL AND PANCREATIC MARKERS	METHOD	REGULATORY
DZ062A	N-Acetyl-β-Glucosaminidase (NAG) Assay	(Three Vial Liquid Stable) (Colorimetric)	☑️ USA: For Research Use Only*
DZ072B	Creatinine Liquid Reagent Assay	(Dual Vial Liquid Stable) (Enzymatic)	☑️ 510(k) Cleared
DZ132A	Lipase Assay	(Dual Vial Liquid Stable) (Methyl-Resorufin)	☑️ 510(k) Exempt
DZ133C	Cystatin C Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared
CAT. #	SEPSIS MARKERS	METHOD	REGULATORY
DZ558A	Procalcitonin (PCT) Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared Health Canada Registered
DZ626A	Serum Amyloid A (SAA) Assay ♦	(Dual Vial Liquid Stable) (Immunoturbidimetric)	USA: For Research Use Only*
CAT. #	VITAMIN MARKERS	METHOD	REGULATORY
DZ888A	EZ Vitamin D Assay	(Dual Vial Liquid Stable) (Immunoturbidimetric)	☑️ 510(k) Cleared

♦ Coming Soon

*USA: For Research Use Only. Not For Use in Diagnostic Procedures.

FEATURED PRODUCTS